


Hematocrit rotor, 24-place



n = 15,000 min⁻¹
max. RCF 21,382

Cat. No. 1650

Cyto rotor, 6-place



n = 4,000 min⁻¹
max. RCF 2,039

Cat. No. (without suspensions) 1626

Cyto rotor, 8-place


























n = 4,000 min⁻¹
max. RCF 2,218



Cat. No. (with suspensions) 1648

For detailed information on the comprehensive Hettich cyto system, please turn to the pages 156-158.

illustrated with
suspensions 1660
and lids 1661
for cyto chambers

illustrated with
suspensions 1680
for angle chambers

capacity in ml	4	5	6	15	1.1-1.4	2.6-3.4	2.7-3	4.5-5	4.9	7.5-10	10	1.6-5	4-7	8-10	15					
Ø x L in mm	10x88	12x75	12x82	17x100	8x66	13x65	11x66	11x92	13x90	15/16x92	15x102	13x75	13x100	16x100/125	17x120					
Cat. No.	-	0553 ²⁾	0501 ²⁾	0518 ²⁾	blood collection/urine tubes										0509					
 rotor Cat. No. 1615																				
																				
Cat. No.	6305	1054-A	-	-	1054-A			-				1054-A	1058	-	-	1647 ²⁹⁾				
boring Ø x L in mm	11.5x67.5	13.5x60	17.7x88		13.5x60			17.7x88				13.5x60	13.5x79	17.7x88	17x104					
tubes per rotor				12												12	12	12	6	6
max. RCF ²⁾	14,006	13,201	16,582		13,201			16,582				13,201	16,582	16,582	16,582	15,455				
radius in mm	87	82	103		82			103				82	103	103	96					
run-up in sec	40																			
run-down in sec, braked	40																			
temperature in °C ¹⁾	-2																			

standard capillaries, heparinised	basic	self-sealing and mylar-coated
Cat. No.	2074	1071
 rotor Cat. No. 1650	lid as evaluation disk incl. 	
	sealing putty	
Cat. No.	2077	-
boring Ø x L in mm	-	
capillaries per rotor	24	
max. RCF ²⁾	21,382	
radius in mm	85	
run-up in sec	13	
run-down in sec, braked	9	
temperature in °C ¹⁾	-1	

²⁹⁾ Set of 6 carriers. For accommodating conical 15 ml tubes with screw cap, remove 6 opposite carriers from rotor No. 1615 and replace them with adapters No. 1647.