

CERTIFICATE OF ASSURANCE



This certificate is your assurance that samples of TipOne®, TipOne® Repeat & TubeOne® products are tested in an independent laboratory and are certified to be free of any detectable RNase, DNase, DNA or Pyrogens.

Sterile products have been sterilised by e-beam irradiation in accordance with EN and ISO Standards by an independent accredited company.

The following limits are observed:

DNA	< 2 x 10 ⁻⁹ mg
RNase	< 1 x 10 ⁻⁹ Kunitz Units
DNase <	1 x 10 ⁻⁶ Kunitz Units
Endotoxin	< 0.001 EU/ml
e-beam irradiation	min. dose 15kGy (1.5 Mrads)

www.starlabgroup.com

